

The Great Tooth Killer: Epidemic of Cracked Teeth, the Science of Strong Teeth

featuring

 BIOCLEAR

Disclosures:
Dr. Clark has financial interest in Bioclear

A scenic photograph of a sunset over a calm body of water. The sky is filled with vibrant orange and yellow clouds, with the sun low on the horizon. A dark silhouette of a forest line is visible in the middle ground, and its reflection is clearly visible in the still water in the foreground.

For a copy of today's presentations,
you have two choices:

Resource Library — Desktop


The screenshot shows the Bioclear Resource Library desktop interface. At the top, there is a dark navigation bar with the phone number 1.855.712.5327 on the left and links for 'Instructions For Use', 'Contact Us', and 'Account' on the right. Below this is a white header with the Bioclear logo, a search bar, and navigation links for 'Products', 'Education', 'Find A Doctor', 'About Us', 'Library', and a shopping cart icon showing '\$0.00 0 items'. The 'Library' link is highlighted with a red box and a red arrow pointing to it. Below the header is the 'Bioclear Resource Library' section, which includes a search bar for resources, a 'Sort by' dropdown, and a row of filter buttons: 'All', 'Articles', 'Free', 'Presentations', 'Tips & Tricks', 'Videos', and 'Webinars'. The 'Presentations' button is highlighted with a red box and a red arrow pointing to it. Below the filters is the 'Presentations' section, which displays three presentation cards. The first card is titled 'Simply Better Composite Restorations' by David Clark DDS. The second card is titled 'The Great Tooth Killer: Epidemic of Cracked Teeth, the Science of Strong Teeth' by David Clark DDS (@bioclearmatrix.com). The third card is titled 'Modern Composite Dentistry: CRUSH your next Class II' by David Clark DDS and Lauren Wilson DMD, dated March 18, 2026. A 'Back To Top' button is visible at the bottom right of the presentation grid.

- 1) Go to [bioclearmatrix.com](https://www.bioclearmatrix.com)
- 2) Click "Library" at the top right
- 3) Click "Presentations"

Resource Library — Mobile



1) Scan the QR code



OR type in
bioclearmatrix.com/resources/

2) Fill out the form

12:06 PM

BIOCLEAR

Home > Resource Library Form

Resource Library Form

FILL OUT THE FORM BELOW TO GET STARTED

After filling out the form you will be redirected to the Resource Library page.

Bioclear Resource Library

Participant's Name *(Required)*


First

Last



Practice Name

AA bioclearmatrix.com

3) Automatically redirects to Resource Library



Dental schools/GPR integrating the Bioclear Method





DC BIOCLEAR LEARNING CENTER

Tacoma USA · Solihull UK
Varberg Sweden · Cairo Egypt
Syracuse Italy · Taubate Brazil
Livermore CA (Bioclear pediatrics)
Seoul Korea · Madrid/Barcelona
Sydney Australia · Baghdad Iraq



The **Bioclear Learning Center** provides a new pathway for dentists and their office to transition from traditional **GV Black dentistry**, to the extraordinary world of **modern direct dentistry**.



THE SCIENCE OF STRONG RESTORED TEETH



DR. DAVID CLARK



DR. ALEX FOK



Dr Alex Fok is a Mechanical Engineer with expertise in solid mechanics, structural analysis and mathematical modeling. He is currently the Academic Director of the Minnesota Dental Research Center for Biomaterials and Biomechanics (MDRCBB). The MDRCBB is an international leader in the development and application of novel characterization techniques for dental biomaterials, with long-standing collaboration with dental materials manufacturers. Dr Fok's research focuses on the development of techniques for material characterization, nondestructive examination, lifetime prediction and shape optimization of dental restorations. A principal aim of his research is to instill more engineering principles and analytical techniques into the design and assessment of dental restorations and treatments so as to improve their longevity and effectiveness.

In this case-based lecture
today we explore:

- Coronal (vertical and cuspal)
Fracturing
- Snap-Off Fracturing

A close-up clinical photograph of a tooth with a snap-off fracture. The fracture is a clean, horizontal break in the enamel and dentin, exposing the underlying pulp chamber. The pulp chamber is filled with a dark, necrotic material, likely a cast resin core. The surrounding gingiva is pink and appears slightly inflamed. The text "Snap-Off Fracturing" is overlaid in a dark grey box with light blue text.

Snap-Off Fracturing

Today's Summary


- Why things break...
- Modern cavity preparations
- Injection Overmolding as a 3rd option
- FEA of loading of teeth & composite
- Long term outcomes & case studies of composite overlays for cracked teeth
- When to endo, when to extract.

Modern Approach to Composite Restorations

Posterior
Restoration



Black Triangle





Glossary

- Bending
Mode of loading that creates curvature, resulting in tension on one side and compression on the other
- Compression
Squeezing of an element
- Hooke's law
This is obeyed when the strain of a material is proportional to the stress applied
- Shear
Force or displacement that results in an angular change in shape of an element
- Spring
Most fundamental structural element in solid mechanics
- Strain
Displacement per unit length
- Strength
Maximum stress required to cause fracture or debonding
- Stress
Force per unit area
- Stress concentrator
A structural feature that results in raised stresses within a small region
- Structural analysis or design triangle
Material, geometry and load
- Tension
Stretching of an element
- Toughness
Energy required to create unit surface area
- Shape Optimization
- Modulus
- Brittle Materials vs. Ductile and Tough Materials
- Surface Resilience
- Cyclic Fatigue

Epidemic of Cracked Teeth



S, Max. Principal
(Avg: 75%)



2
3

ABAQUS Structural Analysis with Job Name = model-cdm-nl.odb Created on Date = 25-A
ODB: model-cdm-nl.odb ABAQUS/STANDARD Version 6.5-1 Wed May 03 13:28:29 GMT Daylight

Step: Step-1
Increment 0: Step Time = 0.000
Primary Var: S, Max. Principal
Deformed Var: U Deformation Scale Factor: +1.000e+01



Traditional Class II: Before



Bioclear Rejuvenation: After

Courtesy Dr. Charlie Regalado

Case Study “Joe”: Contralateral bicuspids fracture nine years apart

2012: Catastrophic fracture tooth #12. Endo
tooth w/ traditional Class II Resin Restoration




2021: Catastrophic fracture tooth #5. Virgin
Tooth



Case Study 1 “Joe”: Contrateral bicuspids fracture nine years apart

2012: Catastrophic
fracture tooth #12.
Endo tooth w/
traditional Class II
Resin Restoration



Case Study Joe bilateral bicuspid fractures



2005: Catastrophic fracture tooth #12. Endo tooth w/ traditional Class II Resin Restoration



Case Study Joe bilateral bicuspid fractures

2021: Catastrophic fracture of #5. Virgin Tooth



Case Study Joe bilateral bicuspids fractures

2021: Catastrophic fracture #12. Virgin Tooth



Case Study “Joe”: Discussion

2012: Catastrophic fracture tooth #12. Endo tooth w/ traditional Class II Resin Restoration



2021: Catastrophic fracture tooth #4. Virgin Tooth



Case Study “Sue”: Traditional Parallel Walled Cavity Preparation (PWCP) 10 years post operative

2021: Asymptomatic tooth at New Patient Examination


When the tooth fractures: Crack initiates at the central groove of composite or at interfacial bond.



When the tooth fractures: Crack initiates through the central groove of composite or at interfacial bond?



When tooth fractures: Crack initiates through the **central groove of composite** or at interfacial bond?



When the tooth fractures: Crack initiates through the central groove of composite or at **interfacial bond**?



When the tooth fractures: Crack initiates through the central groove of composite or at **interfacial bond**?



A close-up photograph of a human tooth showing a significant crack and a grayish discoloration on its surface. The crack runs vertically down the center of the tooth. The surrounding gum tissue is pink and appears slightly inflamed. The text is overlaid in white, bold, sans-serif font.

If the tooth is
gray, the tooth is
cracked

According to the FEA, what two things need to happen to eliminate fracturing?

- 1) No Central Groove
- 2) Horizontalize the prep to create compression (eliminate vertical walls with massive radius bevel in enamel, the prep will look like a trumpet or a calla lily)



Studies Supporting the Bioclear Method

Complete

- Comparing Conventional to Saucer-Shaped Cavity Designs

Dr. Alex Fok, BEng, PhD, MSc

Dr. Hooi Pin Chew, BDS, PhD, FDSRCS

MN Dental Research Center for Biomaterials and Biomechanics

- Comparison of Class II Adaptation and Placement Times

Dr. Richard Price, BDS, DDS, MS, FDS RCS, FRCD(C), PhD

Dept. of Clinical Dental Sciences & Biomedical Engineering Dalhousie University

- Effect of Preheating/Fatiguing/Thermocycling on Mechanical Properties

Taiseer A. Sulaiman, DDS, PhD



Assistant professor, Division Director of Operative Dentistry and Biomaterials, UNC School of Dentistry

- 3M Extraction and Pulp Temperature Testing

Brad Bagley, PhD, DABT Advanced Toxicology Specialist

- 3M Material Property Testing Including Injection Molding

Timothy D. Dunbar, Ph.D. Advanced Product Development Specialist



In Process

- Biofilm Adhesion Study



Sabrina F. Sochacki, DDS, MS, PhD

Indiana University School of Dentistry

Basics of Solid Mechanics





All solid structures can be considered as an assembly of many interconnected springs.






Enamel vs. Dentin Bonding

Bonding to enamel is mostly micromechanical.




Bonding to dentin is more complicated, involving a hybrid layer and exposed collagen that are prone to enzymatic attack.

Anisotropy : tooth tissues




From Orban (1949) Oral Histology and Embriology, St. Louis: C.V. Mosby.



No bending stresses

Isotropy: resin composite



VS.



JOURNAL OF ESTHETIC AND RESTORATIVE DENTISTRY

Official Publication of the American Academy of Esthetic Dentistry,
Japan Academy of Esthetic Dentistry, International Federation of Esthetic Dentistry,
American Academy of Cosmetic and Adhesive Dentistry,
Belgian Academy of Esthetic Dentistry, Dutch Academy of Esthetic Dentistry,
and the Scandinavian Academy of Esthetic Dentistry



www.blackwell.com

Definitive Diagnosis of Early Enamel and Dentinal Cracks Based on Microscopic Evaluation

DAVID J. CLARK, DDS*
GERRYLEE G. SHEETS, DDS†
JACINTHE M. PAQUETTE, DDS‡

ABSTRACT

The diagnoses of cracked teeth and incomplete coronal fractures have historically been symptom based. The dental operating microscope at $\times 16$ magnification can fundamentally change a clinician's ability to diagnose such conditions.

Clinicians have been observing cracks under extreme magnification for nearly a decade. Patterns have become clear that can lead to appropriate treatment prior to symptoms or dissemination to tooth structure occur. Conversely, many cracks are not structural and can lead to misdiagnosis and overtreatment. Methodic microscopic examination, an understanding of crack progression, and an appreciation of the types of cracks will guide a doctor to making appropriate decisions.

Teeth can have structural cracks in various stages. To date, diagnosis and treatment are very often at end stage of crack development.

CLINICAL SIGNIFICANCE

This article gives new guidelines for recognition, visualization, classification, and treatment of cracked teeth based on the routine use of $\times 16$ magnification. The significance of enamel cracks as they relate to dentinal cracks is detailed.

(*J Esthet Restor Dent* 15:XXX-XXX, 2003)

Macroscopic and symptom-driven diagnoses have been the accepted modalities for cracked teeth. The inherent limitation of the lack of visual confirmation result in therapies that often come too late in the treatment process. One lacking first impression of vision through a clinical microscope is the staggering array of cracks that exist within tooth structures. Traditional visualization limited

or ocular astigmaty limits the clinician's ability to assess the presence or severity of the majority of these cracks (Figure 1).

At extreme magnification levels ($\times 14$ and greater), the translucent nature of enamel yields a wealth of information. Subtle color changes within the enamel may indicate early decay, microleakage, and a lack of structural integrity of dentin and

enamel. Being able to see previously invisible crazes can lead restorative dentists to more appropriate early treatment of compromised teeth before devastating fractures, pulpal involvement, and pericoronal breakdown occur. The value of early diagnosis of the structural breakdown of teeth will become even more significant with our aging population coupled with increased tooth retention in this population.

*President, Academy of Microscopic Enhanced Dentistry
†Governator Dentist, Newport Coast Oral Facial Institute, Newport Beach, CA, dental professor,
Restorative Dentistry, USC School of Dentistry, Los Angeles, CA, USA
‡Governator Dentist, Newport Coast Oral Facial Institute Newport Beach, CA, associate professor,
Restorative Dentistry, USC School of Dentistry, Los Angeles, CA, USA

Nomenclature and Classification System for Enamel Cracks

Type I: Have little or no risk of underlying pathology

Type II: Have moderate risk of underlying pathology

Type III: Have high risk of underlying pathology

Type I: Have little or no risk of underlying pathology

- A) Craze lines
- B) Small vertical cracks
- C) Cracks that follow natural anatomic grooves
- D) Cracks with superficial stain penetration

Type II: Have moderate risk of underlying pathology

- E) V-shaped enamel ditching with or without an adjoining restoration often associated with a wear facet centered over an otherwise benign crack
- F) Cracks that detour from or do not follow natural anatomic grooves

Type III: Have high risk of underlying pathology

- G) Diagonal cracks branching off a vertical crack
- H) Horizontal or diagonal cracks that emanate from the corner of a restoration
- I) Cracks that house debris

Type III: Have high risk of underlying pathology

- J. Pair of cracks that outline an area (cusp or marginal ridge) of discolored enamel
- K. Crack with corresponding halo of brown, grey or white centered on crack
- **The tooth or the cusp is gray**


**WHERE ARE WE WITH
CAVITY SHAPES?**

1890

**G. V. BLACK,
M.D., D.D.S.,
Sc.D., LL.D.,
1836-1915**



**MORTICE AND
TENON JOINTS
NECESSARY FOR
MECHANICAL
RETENTION OF
AMALGAM**





EDGETAPER
PLATINUM™



PER GOLD® \$67/6pk
Taper Platinum™

95/6PK

ROTAPER GOLD®
ROTAPER®

FIREWIRE™ NITI

ance as PROTAPER GOLD®
ance as PROTAPER®

THE JOURNAL OF MULTIDISCIPLINARY CARE

Decisions™ IN DENTISTRY

CONSERVATIVE APPROACH TO DEEP CARRIES LESIONS

NATHANIEL LAWSON, DMD, PhD, and
AUGUSTO ROBLES, DMD, MS


RISK FACTORS IN
ORAL CANCER
SCREENING
SARAH GLASS, DDS

HOW MANAGED
CARE WILL IMPACT
DENTAL PRACTICE
JOEL H. BERG, DDS, MS


FEBRUARY 2019

DecisionsInDentistry.com

A Peer-Reviewed Journal






Is this the
year 2026
or 1890?



Stress/Strain Concentration

- Abrupt changes in geometry
- Mismatches in mechanical properties
- Concentrated loads




Epidemic of Cracked Teeth...Why?


**Brittle Materials with
Tension joints**

Composite Resin Material Properties

Compressive vs Diametral Tensile Strength




Compression vs. Tension is “Engineering 101”



compression



tension


A microscopic view of a tooth with a composite filling. The background is a light blue, textured surface representing the tooth's enamel. A dark, rectangular area in the center represents the composite filling. The text is overlaid on a black rectangular background in the center of the image.

But composite
strengthens a
tooth, right?

Does composite strengthen the tooth?


No more than does amalgam

Wahl MJ, Schmitt MM,
Overton DA, Gordon MK.
Prevalence of cusp
fractures in teeth restored
with amalgam and with
resin-based composite.
JADA 135:1127-1132
(2004)










Boxes are bad.

Bevels make them less bad.

The Clark Class II replaces
beveling with radius walls.



It's not the filling, it's the
prep that matters.

The New Science of Strong Materials

or Why You Don't Fall through the Floor

J. E. Gordon



The Epidemic of Cracked and Fracturing Teeth



MAGNIFICATION

Tooth and endodontic failure are on the rise. A recent study by the American Association of Endodontics (AAE) and the American Dental Association (ADA) found that 10% of teeth in the United States are cracked or fractured. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured.

REVERSE CAUSING AND PREVENTION
The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured.

Phase 1	Phase 2	Phase 3	Phase 4
Initial and final diagnosis	Diagnosis of all teeth	Diagnosis of all teeth	Diagnosis of all teeth
Diagnosis of all teeth	Diagnosis of all teeth	Diagnosis of all teeth	Diagnosis of all teeth
Diagnosis of all teeth	Diagnosis of all teeth	Diagnosis of all teeth	Diagnosis of all teeth

- 1) Initial and final diagnosis
- 2) Diagnosis of all teeth
- 3) Diagnosis of all teeth
- 4) Diagnosis of all teeth



Figure 6. An alternative type of crown can be used to restore a cracked tooth. The crown is made of a material that is stronger than the tooth, and the crown is made of a material that is stronger than the tooth.

RESTORATIVE

Fracture Resistant Endodontic and Restorative Preparations



Figure 1. Photograph of the tooth showing the endodontic and restorative procedures. The tooth is shown in various stages of preparation and restoration.



INTRODUCTION
The purpose of this study is to evaluate the effectiveness of various endodontic and restorative preparations. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured.

The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured.

CONCLUSION
The study found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured. The study also found that 10% of teeth are cracked or fractured.

Restoratively, It's Not the Size of the Hole, but the Shape of the Hole

The formula for determining stress concentration due to a cavity preparation is $(1+2\sqrt{\frac{L}{R}})$ where L is length of the cut and R is the radius of the cut. In simple terms, the longer the cut, the worse the cut. A long narrow cut, interestingly, is worse than a wide, round-bottomed cut. One reason that intracoronal composites do not protect the tooth from long-term from


A microscopic image of a metal surface, likely steel, showing a sharp V-shaped notch and a crack extending from its base. The surface has a fine, granular texture. The lighting is bright, highlighting the edges of the notch and the crack.

The formula for stress concentration in a cavity is


$$1 + 2 \sqrt{L/R}$$

L is the depth of the cut.

R is the radius of the cut.

A close-up photograph of a diamond being cut. The diamond is the central focus, showing its facets and the cutting process. The background is a dark, textured surface, likely the cutting table. The lighting is bright, highlighting the facets of the diamond.


Simply: the narrower, the
boxier, and the longer the cut,
the worse the cut.

A microscopic image of a material surface, likely a metal or ceramic, showing a crack. The crack is a dark, irregular line running diagonally across the frame. The surface has a granular texture with various small features and imperfections. The lighting is bright, creating a strong contrast between the dark crack and the lighter surrounding material.


Basic mechanical engineering principle:
The best way to delay crack propagation in a brittle material
is to mitigate a tension area with a compression joint.

And don't unnecessarily connect the dentinal dots


The Prep?



Final views of completed Bioclear Class II Preparation.



To eliminate mechanical retention, we need 2 mm of enamel in compression at the primary strike point of occlusion



Modern Restoration

Design=**I.C.E.**

Infinity edge


Compression based

Enamel driven

Modern Method for Composite Restorations

Clear Anatomic Matrices

- Anterior & Posterior Matrices designed to mimic nature



Injection Mold Composite

- Injection mold warmed Restorative materials
- Industry leading polish, esthetic, strength & wear



Preparation Design

- Designed for composite
- Minimizes stress concentration
- Maximizes enamel involvement



Biofilm Removal

- Remove biofilm before bonding
- Allows bonding to uncut enamel
- Allows infinity edge margins



Systematic restorative protocol for esthetic long-term clinical outcomes

Final Polish

- 3M™ Sof-Lex™ XT coarse discs for reduction
- “Rock Star” polish with Bioclear Magic Mix & RS Polisher



Clark Class II Prep Design



Unique Aspects of the Clark Class II Prep Design




- Removal of biofilm

Unique Aspects of the Clark Class II Prep Design



- Removal of biofilm
- Rounded internal line angles

Unique Aspects of the Clark Class II Prep Design



- Removal of biofilm
- Rounded internal line angles
- Rounded external line angles via the ***Radius Bevel***

Unique Aspects of the Bioclear Prep Design



- Removal of biofilm
- Rounded line angles
- Injection mold rather than place incrementally

Unique Aspects of the Clark Class Prep Design



- Removal of biofilm
- Rounded line angles
- Rounded external line angles via the *Radius Bevel*
- Infinity Edge margins
 - Maximized enamel rod Engagement *with Large Areas of Additive Dentistry*



before

The Bioclear method





**Infinity edge
of the T.R.I.**

The Bioclear Injection Molding Approach



Failing composite;
traditional preparation



Re-restored using the
Bioclear approach



Step by Step Guide for Injection Molded Class II

Studies Supporting the Bioclear Method

Complete

- Comparing Conventional to Saucer-Shaped Cavity Designs

Dr. Alex Fok, BEng, PhD, MSc

Dr. Hooi Pin Chew, BDS, PhD, FDSRCS

MN Dental Research Center for Biomaterials and Biomechanics

- Comparison of Class II Adaptation and Placement Times

Dr. Richard Price, BDS, DDS, MS, FDS RCS, FRCD(C), PhD

Dept. of Clinical Dental Sciences & Biomedical Engineering Dalhousie University

- Effect of Preheating/Fatiguing/Thermocycling on Mechanical Properties

Taiseer A. Sulaiman, DDS, PhD



Assistant professor, Division Director of Operative Dentistry and Biomaterials, UNC School of Dentistry

- 3M Extraction and Pulp Temperature Testing

Brad Bagley, PhD, DABT Advanced Toxicology Specialist

- 3M Material Property Testing Including Injection Molding

Timothy D. Dunbar, Ph.D. Advanced Product Development Specialist



In Process

- Biofilm Adhesion Study

Sabrina F. Sochacki, DDS, MS, PhD

Indiana University School of Dentistry



Fracture Simulation of Class-I Restorations





Prepared by: Ning Ye
Date: 1/11/2019

Reference no. MDRRCBB-NY-JAN11-2019

Traditional Class I Prep



Calla Lily Class I Prep



$$A_1 N_1^B = A_2 N_2^B \frac{N_1}{N_2} = \left(\frac{A_2}{A_1} \right)^{1/B}$$


...if Restoration 1 has a fast fracture strength that is **10% higher than Restoration 2**, then its fatigue life is roughly 3 times as long.

If the difference in fast fracture strength is 20%, then Restoration 1 will last roughly 8 times as long, etc.


Alex Fok

Minnesota Dental Research Center
for Biomaterials and Biomechanics

Calla Lily: moderate occlusal caries or re-treatment of Traditional prep




Traditional Class I filled with composite resin




The compressive strength of composite is about 4x the diametral tensile strength.

Tensile



Compressive



Bioclear Learning Center Guidelines:

- In compression, resin over enamel has no minimum thickness requirement
- In compression, resin over dentin needs at least 2mm thickness
- For Incisal Edges and cusp tips we need 2 mm (occluso gingivally and buccal-lingual)
- All monolithic-no seams or grooves

Calla Lily Prep


Calla Lily Prep






Enamel

Dentin





Enamel

Dentin




Do you trust enamel bonding or **dentin bonding**?

Enamel





Enamel














At what point will this tooth go from compression into tension, and vice versa?



Shape Optimization:
As operators we can *flatten cusps* or *fill in the fossae* (or both).



Enamel vs. resin composite




Enamel rods?


Resin composite



Shape Optimization:
Operators we can ***flatten cusps*** or ***fill in the fossae*** (or both).




Shape Optimization:
Build a bridge from buccal to lingual




Is there anything illegal or immoral
about the anatomy of a maxillary
first molar?


What unique trait does a maxillary first molar possess that protects it from cracking?



The oblique ridge



The oblique bridge



OUCH!!!!!!



Man Copies Nature's Mistake:

Vertical Layering and Deep Continuous Grooves in a Brittle Material (Composite or Ceramic)










5 year post-op fracture of monolithic zirconia

Stress/Strain Concentration

- Abrupt changes in geometry
- Mismatches in mechanical properties
- Concentrated loads



Symptomatic Cracked Lower Molars: Bioclear Consult (2nd Opinion)

Julia's Seattle dentist has treatment
planned her for at least one root
canal and crown **#30 and #31**


Meet Julia











Regardless of symptoms, there is **95%** chance that the first molar is cracked, and only a 20% chance that the 2nd molar is cracked



Regardless of symptoms, there is **95%** chance that the first molar is cracked, and only a **20%** chance that the 2nd molar is cracked





Julia: The wet cotton roll test **benefits:**

- **Trust:** Patient owns the diagnosis
- **Time:** Patients are often confused about which tooth is the problem, this saves time consuming and trust degrading debate
- **Efficiency and Prevention:** You may find two or more cracked teeth

Julia: The wet cotton roll test:

- For reversible pulpitis (RP) cases
- Irreversible pulpitis (IP) gets the percussion test
- Get one end of the cotton roll wet
- Patient holds on to the dry end
- Start at the canine (a negative control)
- Repeat the full test at least once

Julia: The patient driven wet
cotton roll test




Julia: The wet cotton roll test:

- For reversible pulpitis (RP) cases (irreversible pulpitis (IP) gets the percussion test)
- Get one end of the cotton roll wet
- Patient holds on to the dry end
- Start at the canine (a negative control)
- Repeat the full test at least once

Patient Name: Julian R.
Treating Doctor: David Clark
TX Date: 10/23/2023

Lower Arch



OPTION A

Existing restorations will be removed to place the Bioclear restoration. Important step to ensure monolithic color unity and strength to the Bioclear restoration.

- **Disassembly**

- Tooth # 30 \$150
- Tooth # 31 \$150


Posterior Overlay: Replace existing failing restorations, cracked tooth, rejuvenation of tooth appearance and color.

- **Bioclear Posterior Overlay**

- Tooth # 30 \$836
- Tooth # 31 \$836

Patient Name: Julian R.
Treating Doctor: David Clark
TX Date: 10/23/2023

Lower Arch



OPTION A

Existing restorations will be removed to place the Bioclear restoration. Important step to ensure monolithic color unity and strength to the Bioclear restoration.


- **Disassembly**

- Tooth # 30 \$150
- Tooth # 31 \$150

Posterior Overlay: Replace existing failing restorations, cracked tooth, rejuvenation of tooth appearance and color.

- **Bioclear Complex Posterior Overlay**

- Tooth # 30 \$1,200
- Tooth # 31 \$1,200



We show her the patient book and cracked teeth pamphlet. Go straight to page 42. Then she chooses a crown or Bioclear. She will get a little insurance coverage for Bioclear. (We are PPO)

Getting paid to be conservative

BIOCLEAR VS. CROWNS AND VENEERS

Bioclear is an alternative to traditional dental methods for enhancing a smile. Rather than grinding down a tooth to prepare for a crown or veneer, the Bioclear Method of encasing the natural tooth structure retains the natural tooth enamel and tooth durability.



BIOCLEAR PREPARATION

Bioclear allows dentists to **conserve** healthy tooth structure



CROWN PREPARATION

Crowns require dentists to **remove an average of 76%** of the tooth structure prior to the procedure



Julia's overlay treatment images





































Time to *harmonize* the occlusion





Before Harmonizing




After Harmonizing

Remember that composite absorbs 1% of its volume in water in 24 hours, so you need HYPO OCCLUSION on composite



6 week
follow up



2.5 year follow up



Pre-op



What are the **modern occlusal** cavity preps?

1. Fissurotomy
2. Calla Lily
3. Flattened Calla Lily
4. Cuspal Overlay

What are the modern Interproximal cavity preps?

1. Opportunistic Class II
2. Clark Class II...double serpentine,
additive and compression based
3. Clark Class II plus Lateral Partial
Cuspal Overlay













Pre-op



1 year follow up



1 year follow up



2 year follow up



7 year follow up



9 year follow up



9 year follow up



12 year follow up



Pre-op




12 year follow up




An Extreme Calla Lily Clinical Case

Patient heard a loud *snap* when he bit on a seed and now has severe cold and percussion pain





Stoma on Facial




Normal Probing Depth




Normal Probing Depth











18 month follow up









“Build a bridge from buccal to lingual, and the crack(s) becomes dormant.”



Before



Before

“That’s it”



DR. ALEX FOK




“David that’s not going to work, do a crown.”





7 years post-op

Stress/Strain Concentration

- Abrupt changes in geometry
- Mismatches in mechanical properties
- Concentrated loads



	Elastic Modulus (GPa)	Poisson's Ratio	Effective Shrinkage Strain*
Composite	10	0.3	0.34%
Dentin	18	0.2	NA
Enamel	48-80	0.23	NA
Gold	79		
Zirconia	200		



[1] Alex SL Fok, Dent Mater, 2013

- Why things break...
- Modern cavity preparations
- Injection Molding of composite
- FEA of both load and shrinkage
- Long term outcomes & case studies
- The Bioclear Method as a “3rd option”

■ When to do endo, when to extract

- Why things break...
- Modern cavity preparations
- Injection Molding of composite
- FEA of both load and shrinkage
- Long term outcomes & case studies
- The Bioclear Method as a “3rd option”


■ When to do endo, when to extract

JOE

Journal of Endodontics

December 2007 Volume 33, Number 12

Dec
2007



Effect of Environment on Cyclic Fatigue
Page 1433

Clinical page 1405

Association of Pulpitis in a Cohort of 796 Cracked Teeth

Basic page 1430

Initial Evaluation of Ceramicrete as a Potential Root End Filling Material

Case Report page 1484

Nerve Damage Associated with Overextended Obturating Material



Official Journal of the American Association of Endodontists

A Six Year Evaluation of Cracked Teeth Diagnosed with Reversible Pulpitis: Treatment and Prognosis

Keith V. Krell, DDS, MS, MA,* and Eric M. Rivera, DDS, MS†

Abstract

The purpose of this investigation was to report on the clinical outcomes of cracked teeth diagnosed with reversible pulpitis (RP). Eight thousand one hundred seventy-five patients referred for evaluation and treatment during a 6-year period had medical and dental histories, radiographs, pulpal and periapical diagnosis, periodontal probings, direct identification of crack(s) with transillumination, and biting responses on various cusps recorded. All data were stored daily in a database. All cases were treatment planned according to the pulpal and periapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of periapical diagnosis. All patients were recalled at 1 year unless root canal treatment was needed before the anniversary. Results indicated that cracks were identified in 9.7% (796 of 8175) of all teeth evaluated during this time period. Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. To date, none of the original remaining 100 cases of RP have required root canal treatment. The outcomes of this study suggest that if a marginal ridge crack is identified early enough in teeth with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. (*J Endod* 2007;33:1405–1407)

Key Words

Cracked teeth, reversible pulpitis

From the *Department of Endodontics, College of Dentistry, University of Iowa, Iowa City, Iowa; and †Department of Endodontics, School of Dentistry, University of North Carolina, Chapel Hill, North Carolina.

Address requests for reprints to Keith V. Krell, DDS, MS, MA, 1450 28th St, West Des Moines, IA 50266. E-mail address: krekth@iwoi.com.

0099-2595/07 - see front matter

Copyright © 2007 by the American Association of Endodontists.

doi:10.1016/j.joen.2007.08.015

The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" (<http://www.aae.org/NR/rdonlyres/7D73B05C-FEE4-4B00-AB57-086056F163BC/0/bw97ecfc.pdf>) defined 5 types of tooth cracks. Four of the 5 cracks are associated with coronal defects generated from chewing and biting events (1). Teeth with craze lines have no pain, show lines in the enamel, but no "shadows" in the dentin with transillumination. Teeth with fractured cusps have mild pain to biting on a specific cusp, usually a marginal ridge and buccal or lingual groove crack in the dentin, seen as a shadow with transillumination and a Class II restoration. Removal of the restoration might result in the cusp breaking off. Cracked teeth might or might not have a restoration, will exhibit acute pain on mastication, early brief pain to cold, centrally located mesial-to-distal marginal ridge crack(s) seen in the dentin as a shadow(s) with transillumination, normal-to-deep periodontal probings associated with the crack, no detectable movement of cusps with an explorer, and might require removal of existing restorations to definitively diagnose. The pulpal and periapical diagnoses are dependent on the extent of the crack and duration of the symptoms. Split teeth usually have marked pain to chewing, can be considered an extension of the cracked tooth, have mesiodistal cracks extending across both marginal ridges with detectable buccal and lingual cusp separation with an explorer, and deep probings with both marginal ridges. Vertical root fractures begin in the roots of teeth that usually have had endodontic treatment, extend in the buccal-to-lingual plane, usually have minimal signs or symptoms, might have normal probings, and require surgical assessment to complete the diagnosis.

The treatment of teeth diagnosed as a cracked tooth has largely been variable and empirical. The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" suggested "the treatment plan will vary depending on the location and extent of the crack" and noted that "any thermal sensitivity probably indicates the crack extends near or into the pulp, and root canal treatment will be necessary prior to restoring the tooth with a crown" (1). The juxtaposition to assuming the tooth will need root canal treatment before the crown is to place the crown first, see whether symptoms subside, and only perform root canal treatment when the pulpal and periradicular status dictates treatment. Ailor (2) presented a "flow chart" that took into consideration the pulp status at the time of discovery of the crack. He suggested temporizing the tooth with a temporary crown and monitoring it for symptoms.

Incidence data from Hiatt (3) and Cameron (4) found that the most frequently involved teeth were mandibular molars followed by maxillary premolars and then by maxillary first molars. Weine et al (5) found that the most frequently involved teeth were mandibular molars followed by maxillary molars and then by maxillary premolars.

To date there have been no studies that have looked at teeth diagnosed with cracks (cracked teeth) that have a pulpal diagnosis of reversible pulpitis (RP) and are subsequently crowned without initial root canal treatment. The purpose of this investigation was to report on the clinical outcomes of cracked teeth that were initially diagnosed with RP that were crowned during a 6-year period. Specific questions to be answered were the following: (1) what are the locations of these cracked teeth in the mouth? (2) What is the periodontal status of these cracked teeth? (3) Where are the cracks located on the teeth? (4) How many teeth required root canal therapy? (5) What clinical changes occurred in or around the teeth requiring root canal therapy?

Dec 2007



A Six Year Evaluation of Cracked Teeth Diagnosed with Reversible Pulpitis: Treatment and Prognosis

Keith V. Krell, DDS, MS, MA,* and Eric M. Rivera, DDS, MS†

Abstract

The purpose of this investigation was to report on the clinical outcomes of cracked teeth diagnosed with reversible pulpitis (RP). Eight thousand one hundred seventy-five patients referred for evaluation and treatment during a 6-year period had medical and dental histories, radiographs, pulpal and periapical diagnosis, periodontal probings, direct identification of crack(s) with transillumination, and biting responses on various cusps recorded. All data were stored daily in a database. All cases were treatment planned according to the pulpal and periapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of periapical diagnosis. All patients were recalled at 1 year unless root canal treatment was needed before the anniversary. Results indicated that cracks were identified in 9.7% (796 of 8175) of all teeth evaluated during this time period. Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. To date, none of the original remaining 100 cases of RP have required root canal treatment. The outcomes of this study suggest that if a marginal ridge crack is identified early enough in teeth with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. (*J Endod* 2007;33:1405-1407)

Key Words

Cracked teeth, reversible pulpitis

From the *Department of Endodontics, College of Dentistry, University of Iowa, Iowa City, Iowa; and †Department of Endodontics, School of Dentistry, University of North Carolina, Chapel Hill, North Carolina.

Address requests for reprints to Keith V. Krell, DDS, MS, MA, 1450 28th St, West Des Moines, IA 50266. E-mail address: krekth@iowai.com.

0099-2595/07 - see front matter

Copyright © 2007 by the American Association of Endodontists.

doi:10.1016/j.joen.2007.08.015

The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" (<http://www.aae.org/NR/rdonlyres/7D73B05C-FEE4-4B00-AB57-086056F163BC/0/bw97ecfe.pdf>) defined 5 types of tooth cracks. Four of the 5 cracks are associated with coronal defects generated from chewing and biting events (1). Teeth with craze lines have no pain, show lines in the enamel, but no "shadows" in the dentin with transillumination. Teeth with fractured cusps have mild pain to biting on a specific cusp, usually a marginal ridge and buccal or lingual groove crack in the dentin, seen as a shadow with transillumination and a Class II restoration. Removal of the restoration might result in the cusp breaking off. Cracked teeth might or might not have a restoration, will exhibit acute pain on mastication, early brief pain to cold, centrally located mesial-to-distal marginal ridge crack(s) seen in the dentin as a shadow(s) with transillumination, normal-to-deep periodontal probings associated with the crack, no detectable movement of cusps with an explorer, and might require removal of existing restorations to definitively diagnose. The pulpal and periapical diagnoses are dependent on the extent of the crack and duration of the symptoms. Split teeth usually have marked pain to chewing, can be considered an extension of the cracked tooth, have mesiodistal cracks extending across both marginal ridges with detectable buccal and lingual cusp separation with an explorer, and deep probings with both marginal ridges. Vertical root fractures begin in the roots of teeth that usually have had endodontic treatment, extend in the buccal-to-lingual plane, usually have minimal signs or symptoms, might have normal probings, and require surgical assessment to complete the diagnosis.

The treatment of teeth diagnosed as a cracked tooth has largely been variable and empirical. The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" suggested "the treatment plan will vary depending on the location and extent of the crack" and noted that "any thermal sensitivity probably indicates the crack extends near or into the pulp, and root canal treatment will be necessary prior to restoring the tooth with a crown" (1). The juxtaposition to assuming the tooth will need root canal treatment before the crown is to place the crown first, see whether symptoms subside, and only perform root canal treatment when the pulpal and periradicular status dictates treatment. Ailor (2) presented a "flow chart" that took into consideration the pulp status at the time of discovery of the crack. He suggested temporizing the tooth with a temporary crown and monitoring it for symptoms.

Incidence data from Hiatt (3) and Cameron (4) found that the most frequently involved teeth were mandibular molars followed by maxillary premolars and then by maxillary first molars. Weine et al (5) found that the most frequently involved teeth were mandibular molars followed by maxillary molars and then by maxillary premolars.

To date there have been no studies that have looked at teeth diagnosed with cracks (cracked teeth) that have a pulpal diagnosis of reversible pulpitis (RP) and are subsequently crowned without initial root canal treatment. The purpose of this investigation was to report on the clinical outcomes of cracked teeth that were initially diagnosed with RP that were crowned during a 6-year period. Specific questions to be answered were the following: (1) what are the locations of these cracked teeth in the mouth? (2) What is the periodontal status of these cracked teeth? (3) Where are the cracks located on the teeth? (4) How many teeth required root canal therapy? (5) What clinical changes occurred in or around the teeth requiring root canal therapy?

6 year
evaluation
.....the first
study of
its kind



A Six Year Evaluation of Cracked Teeth Diagnosed with Reversible Pulpitis: Treatment and Prognosis

Keith V. Krell, DDS, MS, MA,* and Eric M. Rivera, DDS, MS†

Abstract

The purpose of this investigation was to report on the clinical outcomes of cracked teeth diagnosed with reversible pulpitis (RP). Eight thousand one hundred seventy-five patients referred for evaluation and treatment during a 6-year period had medical and dental histories, radiographs, pulpal and periapical diagnosis, periodontal probings, direct identification of crack(s) with transillumination, and biting responses on various cusps recorded. All data were stored daily in a database. All cases were treatment planned according to the pulpal and periapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of periapical diagnosis. All patients were recalled at 1 year unless root canal treatment was needed before the anniversary. Results indicated that cracks were identified in 9.7% (796 of 8175) of all teeth evaluated during this time period. Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. To date, none of the original remaining 100 cases of RP have required root canal treatment. The outcomes of this study suggest that if a marginal ridge crack is identified early enough in teeth with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. (*J Endod* 2007;33:1405–1407)

Key Words

Cracked teeth, reversible pulpitis

From the *Department of Endodontics, College of Dentistry, University of Iowa, Iowa City, Iowa; and †Department of Endodontics, School of Dentistry, University of North Carolina, Chapel Hill, North Carolina.

Address requests for reprints to Keith V. Krell, DDS, MS, MA, 1450 Hawkins Drive, Denton, TX 76205. E-mail address: krellk@uioa.edu
0099-2398
Copyright © 2007 by Lippincott Williams & Wilkins
DOI: 10.1097/00006583-200712000-00005



The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" (<http://www.aae.org/NR/rdonlyres/7D73B05C-FEE4-4B00-AB57-086056F163BC/0/w97ecfe.pdf>) defined 5 types of tooth cracks. Four of the 5 cracks are associated with coronal defects generated from chewing and biting events (1). Teeth with craze lines have no pain, show lines in the enamel, but no "shadows" in the dentin with transillumination. Teeth with fractured cusps have mild pain to biting on a specific cusp, usually a marginal ridge and buccal or lingual groove crack in the dentin, seen as a shadow with transillumination and a Class II restoration. Removal of the restoration might result in the cusp breaking off. Cracked teeth might or might not have a restoration, will exhibit acute pain on mastication, early brief pain to cold, centrally located mesial-to-distal marginal ridge crack(s) seen in the dentin as a shadow(s) with transillumination, normal-to-deep periodontal probings associated with the crack, no detectable movement of cusps with an explorer, and might require removal of existing restorations to definitively diagnose. The pulpal and periapical diagnoses are dependent on the extent of the crack and duration of the symptoms. Split teeth usually have marked pain to chewing, can be considered an extension of the cracked tooth, have mesiodistal cracks extending across both marginal ridges with detectable buccal and lingual cusp separation with an explorer, and deep probings with both marginal ridges. Vertical root fractures begin in the roots of teeth that usually have had endodontic treatment, extend in the buccal-to-lingual plane, usually have minimal signs or symptoms, might have normal probings, and require surgical assessment to complete the diagnosis.

The treatment of teeth diagnosed as a cracked tooth has largely been variable and empirical. The Fall/Winter 1997 AAE Colleagues for Excellence article entitled "Cracking the Cracked Tooth Code" suggested "the treatment plan will vary depending on the location and extent of the crack" and noted that "any thermal sensitivity probably indicates the crack extends near or into the pulp, and root canal treatment will be necessary prior to restoring the tooth with a crown" (1). The juxtaposition to assuming the tooth will need root canal treatment before the crown is to place the crown first, see whether symptoms subside, and only perform root canal treatment when the pulpal and periradicular status dictates treatment. Ailor (2) presented a "flow chart" that took into consideration the pulp status at the time of discovery of the crack. He suggested temporizing the tooth with a temporary crown and monitoring it for symptoms.

Incidence data from Hiatt (3) and Cameron (4) found that the most frequently involved teeth were mandibular molars followed by maxillary premolars and then by maxillary first molars. Weine et al (5) found that the most frequently involved teeth were mandibular molars followed by maxillary molars and then by maxillary premolars.

To date there have been no studies that have looked at teeth diagnosed with cracks (cracked teeth) that have a pulpal diagnosis of reversible pulpitis (RP) and are subsequently crowned without initial root canal treatment. The purpose of this investigation was to report on the clinical outcomes of cracked teeth that were initially diagnosed with RP that were crowned during a 6-year period. Specific questions to be answered were the following: (1) what are the locations of these cracked teeth in the mouth? (2) What is the periodontal status of these cracked teeth? (3) Where are the cracks located on the teeth? (4) How many teeth required root canal therapy? (5) What clinical changes occurred in or around the teeth requiring root canal therapy?

127 teeth
with R.P.
were
restored
with crowns
and NO
ENDO

TABLE 1. The Number and Percentage of Cracked Teeth by Tooth Position in the Mouth

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	71	167	46	25	243	231	12	1	796
Total	835	1879	754	549	1380	2100	456	222	8175
%	8.50	8.89	6.10	4.55	17.61	11.00	2.63	0.45	9.74

Materials and Methods

There were 8175 patients included in this study from a private endodontic practice population during a 6-year period. The patients were recorded consecutively as they were referred to the endodontist for evaluation and appropriate treatment during the 6-year period. Besides the standard medical history and subjective history, the endodontist was responsible for the diagnosis of all teeth and recorded the following information for all teeth:

- (1) Pulpal response to cold or hot.
- (2) Periapical response to pressure, palpation, and percussion.
- (3) Buccal and lingual periodontal probings were recorded in the mesial and distal interproximal spaces and furca. These interproximal probings would be directed precisely where marginal ridge cracks were identified to indicate the deepest probing of the crack. A total of 6 probing points were recorded for each tooth.
- (4) Identification of a crack(s) with direct transillumination and visualization with and without magnification. The identified crack had to block light transmission and show a definite shadow with both the buccal and lingual coronal light placement. Teeth not exhibiting a shadow were considered to have "crazings" and were not included in this study.
- (5) Responses to biting on various cusps of the diagnosed tooth, with at least 1 cusp exhibiting pain to biting on either a burlew wheel or Tooth Slooth (Professional Results, Inc, Laguna Hills, Calif).

Teeth were diagnosed with RP if (1) there was no history of spontaneous pain; (2) the response to cold went away in less than 3–5 seconds; (3) there was no radiographic pathology.

No teeth were included in this part of the study that could not be confirmed by visualization as having an identifiable crack even if biting sensitivity was present. Restorations were removed only for patients with pulpal diagnoses that required root canal treatment. Teeth diagnosed as cusp fractures, split teeth, and vertical root fractures were also excluded from this study (6). All patients were recalled at 1 year unless root canal treatment was needed before the anniversary, then recall was 1 year later after treatment. The recall therefore extended into a seventh year of the data collection.

All cases were treatment planned according to the pulpal and periapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of periapical diagnosis.

TABLE 2. The Number and Percentage of Cracked Teeth by Tooth Position with RP Eventually Requiring Root Canal Treatment

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	12	33	8	4	29	41	0	0	127
REV. > IP or NEC.	3	9	0	0	8	7	0	0	27
%	25	27	0	0	28	17	0	0	21

REV, reversible pulpitis; IP, irreversible pulpitis; NEC, necrosis.

Results

Of the 8175 cases seen during the 6-year period, 796 cases were diagnosed as cracked teeth (9.7%). Mandibular second molars (243/796, 30%) had the largest incidence followed by mandibular first molars (231/796, 29%) and maxillary first molars (167/796, 21%). All teeth are included in Table 1.

Cases with RP had the following distribution: mandibular first molars (41/127, 32%), maxillary first molars (33/127, 25%), and mandibular second molars (29/127, 23%). All teeth are included in Table 2.

Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. The distribution of teeth requiring root canal therapy was as follows: maxillary first molars (9/27, 33%), mandibular second molars (8/27, 29%), mandibular first molars (7/27, 26%), and maxillary second molars (3/27, 11%).

All teeth had initial interproximal probings less than 3 mm in the space associated with the identified crack. Increased interproximal probings were associated with the fractured marginal ridge for only 5 of the 27 teeth requiring root canal treatment. The greatest increase in probing depth was 2 mm for 2 of 5 teeth.

The teeth requiring root canal treatment had the crack located on the distal marginal ridge in 15 of the 27 cases (56%). The distribution was mandibular second molars (6/27, 27%), mandibular first molars (3/27, 11%), maxillary first molars (3/27, 11%), maxillary second premolars (2/27, 7%), and maxillary second molars (1/27, 4%). The crack was located only on the mesial marginal ridge in 4 of the 27 cases (15%). The distribution was the maxillary first molars (3/27, 11%) and maxillary first premolars (1/27, 4%). Both marginal ridges were involved in 8 of the 27 cases (29%). The distribution was mandibular first molars (5/27, 19%), maxillary first molars (2/27, 7%), and maxillary second molars (1/27, 4%). None of the teeth had fractures that extended into the floor of the chamber or rendered them "non-restorable".

None of the original remaining 100 cases of RP required root canal treatment.

Discussion

The patients who composed this database were all patients referred to a private practice endodontist. Our incidence data are in agreement with the findings of Weine et al (5), which were also derived from an endodontist's practice. The difference in percentages from

What are the 5 questions you ask to distinguish R.P. ?



3) Has the pain ever woken
you?

4) Lingering pain to heat?

Table 1 Comparison of the performance of the different methods for the detection of the different types of defects.

Defect	Visual	Ultrasonic	Radiographic	Thermal	Acoustic Emission	Strain	Acoustic	Electromagnetic	Other
Cracks	Low	High	High	Low	High	High	High	High	High
Delamination	High	Low	High	High	Low	Low	Low	Low	Low
Porosity	Low	Low	High	Low	Low	Low	Low	Low	Low
Disbond	High	Low	High	High	Low	Low	Low	Low	Low
Corrosion	High	Low	High	High	Low	Low	Low	Low	Low
Impact	High	Low	High	High	Low	Low	Low	Low	Low
Spall	High	Low	High	High	Low	Low	Low	Low	Low
Reinforcement	High	Low	High	High	Low	Low	Low	Low	Low
Delamination	High	Low	High	High	Low	Low	Low	Low	Low
Cracks	Low	High	High	Low	High	High	High	High	High
Delamination	High	Low	High	High	Low	Low	Low	Low	Low
Porosity	Low	Low	High	Low	Low	Low	Low	Low	Low
Disbond	High	Low	High	High	Low	Low	Low	Low	Low
Corrosion	High	Low	High	High	Low	Low	Low	Low	Low
Impact	High	Low	High	High	Low	Low	Low	Low	Low
Spall	High	Low	High	High	Low	Low	Low	Low	Low
Reinforcement	High	Low	High	High	Low	Low	Low	Low	Low
Delamination	High	Low	High	High	Low	Low	Low	Low	Low



5) Is there “radiographic pathosis”? (Krell’s lexicon)

Table 1

Study	Year	Sample Size	Study Design	Outcome
1	1998	100	Case Report	Spontaneous Pain
2	2001	50	Cross-sectional	Pain to heat
3	2003	150	Case Series	Pain at night
4	2005	200	Case Report	Pain to cold
5	2007	300	Cross-sectional	Pain to cold
6	2009	400	Cross-sectional	Pain to cold
7	2011	500	Cross-sectional	Pain to cold
8	2013	600	Cross-sectional	Pain to cold
9	2015	700	Cross-sectional	Pain to cold
10	2017	800	Cross-sectional	Pain to cold



- Spontaneous Pain? Do the endo
- Pain to heat that lingers? Do the Endo
- Pain that wakes you at night? Do the endo
- Pain to cold and chewing that lingers for more than a minute or two? Do the endo

TABLE 1. The Number and Percentage of Cracked Teeth by Tooth Position in the Mouth

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	71	167	46	25	243	231	12	1	796
Total	835	1879	754	549	1380	2100	456	222	8175
%	8.50	8.89	6.10	4.55	17.61	11.00	2.63	0.45	9.74

Materials and Methods

There were 8175 patients included in this study from a private endodontic practice population during a 6-year period. The patients were recorded consecutively as they were referred to the endodontist for evaluation and appropriate treatment during the 6-year period. Besides the standard medical history and subjective history, the endodontist was responsible for the diagnosis of all teeth and recorded the following information for all teeth:

- (1) Pulpal response to cold or hot.
- (2) Periapical response to pressure, palpation, and percussion.
- (3) Buccal and lingual periodontal probings were recorded in the mesial and distal interproximal spaces and furca. These interproximal probings would be directed precisely where marginal ridge cracks were identified to indicate the deepest probing of the crack. A total of 6 probing points were recorded for each tooth.
- (4) Identification of a crack(s) with direct transillumination and visualization with and without magnification. The identified crack had to block light transmission and show a definite shadow with both the buccal and lingual coronal light placement. Teeth not exhibiting a shadow were considered to have "crazings" and were not included in this study.
- (5) Responses to biting on various cusps of the diagnosed tooth, with at least 1 cusp exhibiting pain to biting on either a burlew wheel or Tooth Slooth (Professional Results, Inc, Laguna Niguel, CA).

Teeth were diagnosed with RP if (1) there was no history of spontaneous pain; (2) the response to cold went away in less than 3–5 seconds; (3) there was no radiographic pathology.

No teeth were included in this part of the study that could not be confirmed by visualization as having an identifiable crack even if biting sensitivity was present. Restorations were removed only for patients with pulpal diagnoses that required root canal treatment. Teeth diagnosed as cusp fractures, split teeth, and vertical root fractures were also excluded from this study (6). All patients were recalled at 1 year unless root canal treatment was needed before the anniversary, then recall was 1 year later after treatment. The recall therefore extended into a seventh year of the data collection.

All cases were treatment planned according to the pulpal and periapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of periapical diagnosis.

TABLE 2. The Number and Percentage of Cracked Teeth by Tooth Position with RP Eventually Requiring Root Canal Treatment

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	12	33	8	4	29	41	0	0	127
REV. >IP or NEC.	3	9	0	0	8	7	0	0	27
%	25	27	0	0	28	17	0	0	21

REV, reversible pulpitis; IP, irreversible pulpitis; NEC, necrosis.

Results

Of the 8175 cases seen during the 6-year period, 796 cases were diagnosed as cracked teeth (9.7%). Mandibular second molars (243/796, 30%) had the largest incidence followed by mandibular first molars (231/796, 29%) and maxillary first molars (167/796, 21%). All teeth are included in Table 1.

Cases with RP had the following distribution: mandibular first molars (41/127, 32%), maxillary first molars (33/127, 25%), and mandibular second molars (29/127, 23%). All teeth are included in Table 2.

Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. The distribution of teeth requiring root canal therapy was as follows: maxillary first molars (9/27, 33%), mandibular second molars (8/27, 29%), mandibular first molars (7/27, 26%), and maxillary second molars (3/27, 11%).

All teeth had initial interproximal probings less than 3 mm in the space associated with the identified crack. Increased interproximal probings were associated with the fractured marginal ridge for only 5 of the 27 teeth requiring root canal treatment. The greatest increase in probing depth was 2 mm for 2 of 5 teeth.

The teeth requiring root canal treatment had the crack located on the distal marginal ridge in 15 of the 27 cases (56%). The distribution was mandibular second molars (6/27, 27%), mandibular first molars (3/27, 11%), maxillary first molars (3/27, 11%), maxillary second premolars (2/27, 7%), and maxillary second molars (1/27, 4%). The crack was located only on the mesial marginal ridge in 4 of the 27 cases (15%). The distribution was the maxillary first molars (3/27, 11%) and maxillary first premolars (1/27, 4%). Both marginal ridges were involved in 8 of the 27 cases (29%). The distribution was mandibular first molars (5/27, 19%), maxillary first molars (2/27, 7%), and maxillary second molars (1/27, 4%). None of the teeth had fractures that extended into the floor of the chamber or rendered them "non-restorable".

None of the original remaining 100 cases of RP required root canal treatment.

Discussion

The patients who composed this database were all patients referred to a private practice endodontist. Our incidence data are in agreement with the findings of Weine et al (5), which were also derived from an endodontist's practice. The difference in percentages from

How many pulps eventually died?



TABLE 1. The Number and Percentage of Cracked Teeth by Tooth Position in the Mouth

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	71	167	46	25	243	231	12	1	796
Total	835	1879	754	549	1380	2100	456	222	8175
%	8.50	8.89	6.10	4.55	17.61	11.00	2.63	0.45	9.74

Materials and Methods

There were 8175 patients included in this study from a private endodontic practice population during a 6-year period. The patients were recorded consecutively as they were referred to the endodontist for evaluation and appropriate treatment during the 6-year period. Besides the standard medical history and subjective history, the endodontist was responsible for the diagnosis of all teeth and recorded the following information for all teeth:

- (1) Pulpal response to cold or hot.
- (2) Periapical response to pressure, palpation, and percussion.
- (3) Buccal and lingual periodontal probing were recorded in the mesial and distal interproximal spaces and furca. These interproximal probing would be directed precisely where marginal ridge cracks were identified to indicate the deepest probing of the crack. A total of 6 probing points were recorded for each tooth.
- (4) Identification of a crack(s) with direct transillumination and visualization with and without magnification. The identified crack had to block light transmission and show a definite shadow with both the buccal and lingual coronal light placement. Teeth not exhibiting a shadow were considered to have "crazings" and were not included in this study.
- (5) Responses to biting on various cusps of the diagnosed tooth, with at least 1 cusp exhibiting pain to biting on either a burlew wheel or Tooth Slooth (Professional Results, Inc, Laguna Niguel, CA).

Teeth were diagnosed with RP if (1) there was no history of spontaneous pain; (2) the response to cold went away in less than 3–5 seconds; (3) there was no radiographic pathology.

No teeth were included in this part of the study that could not be confirmed by visualization as having an identifiable crack even if biting sensitivity was present. Restorations were removed only for patients with pulpal diagnoses that required root canal treatment. Teeth diagnosed as cusp fractures, split teeth, and vertical root fractures were also excluded from this study (6). All patients were recalled at 1 year unless root canal treatment was needed before the anniversary, then recall was 1 year later after treatment. The recall therefore extended into a seventh year of the data collection.

All cases were treatment planned according to the pulpal and peritapical diagnosis. Cases with RP were treatment planned for crowns only, regardless of peritapical diagnosis.

TABLE 2. The Number and Percentage of Cracked Teeth by Tooth Position with RP Eventually Requiring Root Canal Treatment

	Maxillary Second Molar	Maxillary First Molar	Maxillary Second Premolar	Maxillary First Premolar	Mandibular Second Molar	Mandibular First Molar	Mandibular Second Premolar	Mandibular First Premolar	Totals
Cracked	12	33	8	4	29	41	0	0	127
REV. >IP or NEC.	3	9	0	0	8	7	0	0	27
%	25	27	0	0	28	17	0	0	21

REV, reversible pulpitis; IP, irreversible pulpitis; NEC, necrosis.

Results

Of the 8175 cases seen during the 6-year period, 796 cases were diagnosed as cracked teeth (9.7%). Mandibular second molars (243/796, 30%) had the largest incidence followed by mandibular first molars (231/796, 29%) and maxillary first molars (167/796, 21%). All teeth are included in Table 1.

Cases with RP had the following distribution: mandibular first molars (41/127, 32%), maxillary first molars (33/127, 25%), and mandibular second molars (29/127, 23%). All teeth are included in Table 2.

Of 127 patients specifically diagnosed with RP, 27 converted to irreversible pulpitis (N = 21) in 58 days or to necrotic pulp (N = 6) in 149 days. The distribution of teeth requiring root canal therapy was as follows: maxillary first molars (9/27, 33%), mandibular second molars (8/27, 29%), mandibular first molars (7/27, 26%), and maxillary second molars (3/27, 11%).

All teeth had initial interproximal probing less than 3 mm in the space associated with the identified crack. Increased interproximal probing were associated with the fractured marginal ridge for only 5 of the 27 teeth requiring root canal treatment. The greatest increase in probing depth was 2 mm for 2 of 5 teeth.

The teeth requiring root canal treatment had the crack located on the distal marginal ridge in 15 of the 27 cases (56%). The distribution was mandibular second molars (6/27, 27%), mandibular first molars (3/27, 11%), maxillary first molars (3/27, 11%), maxillary second premolars (2/27, 7%), and maxillary second molars (1/27, 4%). The crack was located only on the mesial marginal ridge in 4 of the 27 cases (15%). The distribution was the maxillary first molars (3/27, 11%) and maxillary first premolars (1/27, 4%). Both marginal ridges were involved in 8 of the 27 cases (29%). The distribution was mandibular first molars (5/27, 19%), maxillary first molars (2/27, 7%), and maxillary second molars (1/27, 4%). None of the teeth had fractures that extended into the floor of the chamber or rendered them "non-restorable".

None of the original remaining 100 cases of RP required root canal treatment.

Discussion

The patients who composed this database were all patients referred to a private practice endodontist. Our incidence data are in agreement with the findings of Weine et al (5), which were also derived from an endodontist's practice. The difference in percentages from

How many pulps eventually died? 20%



other studies might be largely due to the nature of general dental populations versus an endodontic specialty population. Hiat (5) and Cameron (4) reported maxillary premolars having the second highest incidence of cracks, which represents a periodontist and general dental practice. Therefore, this study's population might be under-represented with respect to premolars and the actual incidence of cracked teeth in the population at large.

This study found that 21% of the cases diagnosed with RP and a crack eventually required root canal treatment. Although no other studies have reported this finding, there have been studies examining the necessity for root canal treatment on teeth restored with crowns. Saunders and Saunders (7) reported that 19% of crowned teeth in a Scottish dental school population had root canal treatment after crown placement. This study was unable to report whether cracks were present before the crown was placed.

Cheung et al (8) found that 15% (19/122) of teeth restored with a ceramo-metallic crown required root canal treatment after crown placement in a population in Asia. They also found that of those serving as an abutment of a fixed-fixed bridge (25/77), 32% required root canal treatment after final cementation of the bridge.

Although the populations in these 2 studies are quite diverse, the similar percentages for teeth requiring root canal treatment after crown placement suggest that 15%–19% of their patients are to be expected for all crowns. The 21% of the cases in this study with RP and cracks that were crowned and subsequently required root canal treatment is only slightly higher than the 15%–19% and suggests similar incidence data.

When examining a tooth with both mesial and distal marginal ridge fractures, the natural assumption would be that more of these teeth

would eventually require root canal treatment. Our data showed more teeth with a single marginal ridge crack, either mesial or distal, eventually required root canal treatment. This underlines the difficulty in predicting the eventual need for root canal treatment in teeth with RP and a cracked marginal ridge.

Conclusions

The outcomes of this study suggest that if a crack is identified early enough in cases with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. Progression of interproximal periodontal defects associated with the crack(s) will occur in a very small percentage of the cases (5/127, 4%).

References

1. Cracking the cracked tooth code. *Endodontics: Colleagues for Excellence* 1997 (Fall/Winter):1–15.
2. Aizer JJ Jr. Margging incomplete tooth fractures. *J Am Dent Assoc* 2000;131:1168–74.
3. Hiat NH. Incomplete crown root fracture in pulpal periodontal disease. *J Periodontol* 1973;44:369–79.
4. Cameron GE. The cracked tooth syndrome: additional findings. *J Am Dent Assoc* 1976;93:971–5.
5. Wene FS, Detharry Jr, James A. Cracked tooth syndrome: vertical fractures of posterior teeth. In: Wene FS, ed. *Endodontic therapy*, 3rd ed. St Louis: Mosby, 1982:8–15.
6. Rivera EM, Walton BE. Longitudinal fractures. In: Torabinejad M, Walton BE, eds. *Principles and practice of endodontics*, 4th ed. Philadelphia: Saunders, in press.
7. Saunders WP, Saunders EM. Prevalence of periradicular periodontitis associated with crowned teeth in an adult Scottish subpopulation. *Br Dent J* 1998;185:137–40.
8. Cheung GS, Lee SC, Ng RP. Fate of vital pulps beneath a metal-ceramic crown or a bridge retainer. *Int Endod J* 2005;38:521–30.

What is the
pulpal death
rate for any
tooth w/ a
crown? 15-
19%



other studies might be largely due to the nature of general dental populations versus an endodontic specialty population. Blatt (3) and Cameron (4) reported mandibular premolars having the second highest incidence of cracks, which represents a periodontist and general dental practice. Therefore, this study's population might be under-represented with respect to premolars and the actual incidence of cracked teeth in the population at large.

This study found that 21% of the cases diagnosed with RP and a crack eventually required root canal treatment. Although no other studies have reported this finding, there have been studies examining the necessity for root canal treatment on teeth restored with crowns. Saunders and Saunders (7) reported that 19% of crowned teeth in a Scottish dental school population had root canal treatment after crown placement. This study was unable to report whether cracks were present before the crown was placed.

Cheong et al (8) found that 15% (19/122) of teeth restored with a ceramo-metallic crown required root canal treatment after crown placement in a population in Asia. They also found that of those serving as an abutment of a fixed-fixed bridge (25/77), 32% required root canal treatment after final cementation of the bridge.

Although the populations in these 2 studies are quite diverse, the similar percentages for teeth requiring root canal treatment after crown placement suggest that 15%–19% of their patients are to be expected for all crowns. The 21% of the cases in this study with RP and cracks that were crowned and subsequently required root canal treatment is only slightly higher than the 15%–19% and suggests similar incidence data.

When examining a tooth with both mesial and distal marginal ridge fractures, the natural assumption would be that more of these teeth

would eventually require root canal treatment. Our data showed more teeth with a single marginal ridge crack, either mesial or distal, eventually required root canal treatment. This underlines the difficulty in predicting the eventual need for root canal treatment in teeth with RP and a cracked marginal ridge.

Conclusions

The outcomes of this study suggest that if a crack is identified early enough in cases with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. Progression of interproximal periodontal defects associated with the crack(s) will occur in a very small percentage of the cases (5/127, 4%).

References

1. Cracking the cracked tooth code. Endodontics: Guidelines for Excellence 1997 (Fall/Winter):1–13.
2. Aher JE Jr. Managing incomplete tooth fractures. J Am Dent Assoc 2000;131:1168–74.
3. Blatt WB. Incomplete crown construction in pulpal periodontal disease. J Periodontol 1973;44:269–70.
4. Cameron CE. The cracked tooth syndrome: additional findings. J Am Dent Assoc 1976;103:971–5.
5. Wetze FS, Desberry Jr, James A. Cracked-tooth syndrome: vertical fractures of posterior teeth. In: Wetze FS, ed. Endodontic therapy. 5th ed. St Louis: Mosby, 1982:6–15.
6. Rivera EM, Wilson RE. Longitudinal fractures. In: Torabinejad M, Walton RE, eds. Principles and practice of endodontics, 4th ed. Philadelphia: Saunders, in press.
7. Saunders WP, Saunders EM. Prevalence of periodontal periodontitis associated with crowned teeth in an adult Scottish subpopulation. Br Dent J 1998;185:137–40.
8. Cheong GS, Lee BC, Ng RP. Fate of vital pulps beneath a metal-ceramic crown on a bridge abutment. Int Endod J 2005;38:521–30.

What is the pulpal death rate for bridge abutments? 32%



other studies might be largely due to the nature of general dental populations versus an endodontic specialty population. Iliatt (3) and Cameron (4) reported mandibular premolars having the second highest incidence of cracks, which represents a periodontist and general dental practice. Therefore, this study's population might be under-represented with respect to premolars and the actual incidence of cracked teeth in the population at large.

This study found that 21% of the cases diagnosed with RP and a crack eventually required root canal treatment. Although no other studies have reported this finding, there have been studies examining the necessity for root canal treatment on teeth restored with crowns. Saunders and Saunders (7) reported that 19% of crowned teeth in a Scottish dental school population had root canal treatment after crown placement. This study was unable to report whether cracks were present before the crown was placed.

Cheong et al (8) found that 15% (19/122) of teeth restored with a ceramo-metallic crown required root canal treatment after crown placement in a population in Asia. They also found that of those serving as an abutment of a fixed-fixed bridge (25/77), 32% required root canal treatment after final cementation of the bridge.

Although the populations in these 2 studies are quite diverse, the similar percentages for teeth requiring root canal treatment after crown placement suggest that 15%–19% of their patients are to be expected for all crowns. The 21% of the cases in this study with RP and cracks that were crowned and subsequently required root canal treatment is only slightly higher than the 15%–19% and suggests similar incidence data.

When examining a tooth with both mesial and distal marginal ridge fractures, the natural assumption would be that more of these teeth

would eventually require root canal treatment. Our data showed more teeth with a single marginal ridge crack, either mesial or distal, eventually required root canal treatment. This underlines the difficulty in predicting the eventual need for root canal treatment in teeth with RP and a cracked marginal ridge.

Conclusions

The outcomes of this study suggest that if a crack is identified early enough in cases with a diagnosis of RP and a crown is placed, root canal treatment will be necessary in about 20% of these cases within a 6-month period. Progression of interproximal periodontal defects associated with the crack(s) will occur in a very small percentage of the cases (5/127, 4%).

References

1. Cracking the cracked tooth code. Endodontics: Colleague for Excellence 1997 (Fall/Winter):1–13.
2. Alier JE Jr. Managing incomplete tooth fractures. J Am Dent Assoc 2000;131:1168–74.
3. Iliatt WB. Incomplete crown construction in pulpal periodontal disease. J Periodontol 1973;44:269–70.
4. Cameron CE. The cracked tooth syndrome: additional findings. J Am Dent Assoc 1970;63:971–5.
5. Wetse FS, Desberry Jr, James A. Cracked-tooth syndrome: vertical fractures of posterior teeth. In: Wetse FS, ed. Endodontic therapy. 5th ed. St Louis: Mosby, 1982:6–15.
6. Rivera EM, Wilson RE. Longitudinal fractures. In: Torabinejad M, Walton RE, eds. Principles and practice of endodontics, 4th ed. Philadelphia: Saunders, in press.
7. Saunders GP, Saunders EM. Prevalence of periodontal periodontitis associated with crowned teeth in an adult Scottish subpopulation. Br Dent J 1998;185:137–40.
8. Cheong GS, Lee BC, Ng RP. Rate of vital pulps beneath a metal-ceramic crown on a bridge retainer. Int J Endod J 2005;38:521–30.

Big picture:
Cracked teeth
with R.P. don't
necessarily
need endo.



Big picture:

Bioclear

Composite

Overlays have
shown immediate
elimination of
symptoms*



Big picture: If you *give up* and do a crown, prepare, impress and temporize normally but wait 6-8 weeks to seat crown



Big picture:

PROBE!

PROBE!

PROBE!



What is our most
reliable clinical test for
irreversible pulpitis?


**PERCUSSION TEST,
and REPEAT**



A Modern Method for Composite Restorations

Anatomic HD Clear Matrices

- Anterior & Posterior Matrices designed to mimic nature



Injection Mold Composite

- Injection mold warmed **Kuraray Majesty anterior composites** and **3M bulk fill** in Class I and Class II restorations
- Like the inverse of my mullet haircut: Party in the front Business in the back

Preparation Design

- Designed for composite
- Minimizes stress concentration
- Maximizes enamel involvement

Bioclear Blaster or Equivalent for Biofilm Removal

- Uses aluminum trihydroxide w/ water
- Allows bonding to uncut enamel
- Allows infinity edge margins



Systematic restorative protocol for esthetic long-term clinical outcomes

Rock Star Polish

- Coarse discs for reduction
- “Rock Star” polish with **Bioclear Magic Mix** and disposable cup
- **Bioclear RS Diamond Polisher**





Step by Step Guide for Injection Molded Class II

Which tooth is cracked?



Which tooth is cracked?



Which tooth is cracked?




Which tooth is cracked?



?






?





**Wear facet so deep
that there is a
marginal ridge
height discrepancy**






Is this pocket
from chronic
food impaction
or from the
crack?



































Learning Objectives from the Case:

- Don't sweep the contacts*
- 3 mm thickness over the crack
- Build a bridge from buccal to lingual
- The crack becomes dormant
- "Harmonize" the opposing w/ consent
- If food impaction continues, so will pain and pocketing
- Contact technique works for other problems**
- We need a powerful separator



Pre-op, post op, and 3 year follow up

∞ BIOCLEAR



Are all
Separators
Equal?

NO!



Stretched to 20mm ONCE!!!





And now most of
the power is
permanently lost

Once the yield of the metal/shape is reached, the
metal undergoes plastic deformation



Stretched to 20mm Twenty Times



Stretched to 20mm Twenty Times




TwinRing Universal






How do we win the *Snug Contact Game?*

- Power
- Strategy
- Technique



How do we win the *Snug Contact Game?*

- Power
- Strategy
- Technique



Nate Lawson DMD PhD


Director of the Division of Biomaterials and Pranit Bora. BDS, MDS.

Resident, Division of Biomaterials UAB School of Dentistry



TwinRing Universal
28 N

TwinRing Universal Instron Comparison



Note: Rings tested are not new and some rings stretch out quickly and lose up to half of their power after multiple uses. 1N = 1kg (m/s/s)

The slip-off test



 BIOCLEAR

Bioclear Matrices - Posterior





EVOLVE --- MATRIX

Premolar



Molar




EVOLVE
MATRIX

In this case based lecture today we explore:

- Coronal (vertical and cuspal) Fracturing
- Snap-Off Fracturing






What Would You Do?





A different
toothache
patient we saw
recently...
























Basics of Solid Mechanics





All solid structures can be considered as an assembly of many interconnected springs.



Stress/strain concentration

This occurs at

- Abrupt changes in geometry
- Mismatches in mechanical properties
- Concentrated loads



What is the percent of
coronal volume removed in a
conservative full porcelain
crown preparation?

~3/4 Crown



53.0

34 *

All Ceramic
.8 Margin



70.4

All Ceramic
1 mm Margin



74.1



All Ceramic
1.4 mm Margin



75.9

Edelhoff D, Sorensen JA.

1 Structure Removal Associated with Various Preparation Designs for Posterior Teeth.
Int J Periodontics Restorative Dent 2002; vol 22; 3: 240-249.



75.9% volume gone + stress concentration

Is monolithic injection molded composite a viable alternative to full ceramic crowns in some cases?

You be the judge of that



Recommended Bioclear Matrices by Indication

For more information contact us
Bioclear Matrix Systems
1-855-712-5327



360 Veneer

Class V and to significantly increase overjet or correct anterior open bite.



TSS Kit

Anterior Teeth:
Class III, Class IV, Class V
Fractured Incisors
Severe Wear
Composite Veneer
Full Composite Crown



BT (Black Triangle) Kit

Anterior Teeth:
Black Triangles
Peg Laterals
Diastema Closure
Instant Ortho
Class V



Evolve Matrix Kit

Posterior Teeth:
Class I, Class II, Class V



TOOTH & SURFACE SPECIFIC MATRIX SYSTEM

BY ∞ BIOCLEAR



140 Anterior Matrices

#6 through #11 in Mesial & Distal
Small & Medium Lower Incisor

75 Wedges

25 Small Wedges
50 Medium Wedges

The next generation of anterior matrices:

We're taking the
guesswork out of to
matrix selection




Contact your Bioclear Sales Rep to Pre-Order

WWW.BIOCLEARMATRIX.COM

1.855.712.5327

The TSS system is for small and large anterior restorations, broken teeth and anterior esthetic veneering when you DO NOT have a black triangle or diastema







Is monolithic injection molded composite a viable alternative to full ceramic crowns in some cases?

You be the judge of that



The patient is a 35- year-old male. His original chief complaint was a discolored filling on the distal of tooth #8. The patient was given two treatment plans, one to simply replace a few defective restorations and remove caries with traditional fillings or in patient terms we said, “We can patch the holes, or I can rejuvenate your smile. The patchwork plan will be healthy but will not make a significant esthetic change. In addition, the severe wear and acid erosion present on the palatal surfaces could eventually lead to catastrophic problems later i.e. root canals, infection, and tooth loss.” The patient opted for Bioclear rejuvenation versus simple fillings because he wanted a beautiful smile. He chose Bioclear in lieu of crowns because he understood that Bioclear is a more conservative and healthier approach to achieving his goals than traditional crowns.

Phase one of the treatment plan was to restore the anterior six teeth, simultaneously opening the vertical dimension to reduce the need for aggressive tooth reduction for material thickness. In addition, because the patient needed to have the teeth lengthened by 2 mm, opening the vertical dimension by 2 mm allowed the overbite to be more ideal. Because the patient could only commit to the cost of the six Bioclear restorations (\$11,600) we placed transitional occlusal flowable composites (thick sealants) on the four maxillary premolar teeth, and we will allow the molars to settle into occlusion utilizing the well-researched Dahl Technique.


Phase two will be to restore the remaining teeth with Bioclear when the patient has his finances ready.

Treatment Summary Short Version:

Treatment was finished in a single 4-hour session. (Less experienced clinicians should plan to give one hour per tooth). Bioclear TSS Matrices were utilized on all the teeth except tooth #10. #10 required the Bioclear BT matrix system to create “instant ortho” and because a diastema was present there. Bioclear Diamond wedges were used as needed in areas where the contact was lost during caries removal or removal of old composites. Bioclear RSP X-course discs (Black) were used to shape the incisal edges and smooth the small seams present where the matrices meet on the facial and palatal. Final polish was achieved with Bioclear Magic Mix and then Rock Star Polish cups and cones.

One-week postoperative visit revealed healthy teeth and gingiva. The patient was ecstatic about his new smile, had zero post-operative pain or sensitivity, and expressed that his new bite with the increase in VDO felt more comfortable than before.






What percentage of the tooth is removed for a conservative crown prep?



Source: Google Images

What percentage of the tooth is removed for a conservative crown prep?



Source: Google Images



Wouldn't it be nice to preserve nearly all the healthy tooth structure and at the same time completely rejuvenate this guy's smile?




Simplified Smile Design!

How long are his centrals currently? 8.5 mm



How wide are his centrals? 8.5 mm



5
10
15
20

Central Width
7.0mm

7.5mm

8.0mm


8.5mm

9.0mm

9.5mm

10mm

 **Panadent**
Esthetic Proportion
Gauge



Central Width

Central Height

8.0mm

10.0mm

8.5mm

10.6mm

9.0mm

11.2mm

9.5mm

11.8mm

10.0mm

12.5mm













TSS

TOOTH & SURFACE SPECIFIC
MATRIX SYSTEM
BY BIOCLEAR



140 Anterior Matrices
45 Preformed Plastic Matrix & Crown
Small & Medium Lower Incisor

75 Wedges
25 Small, Medium
50 Medium Wedges

**The next generation
of anterior matrices:**

We're taking the
guesswork out of to
matrix selection



Contact your Bioclear Sales Rep to Pre-Order
WWW.BIOCLEARMATRIX.COM
1.855.712.5427



TSS

TOOTH & SURFACE SPECIFIC
MATRIX SYSTEM
BY BIOCLEAR



141 Anterior Matrices
10 Standard 1/2 Size Standard 2/3 Size
1/4 Standard 1/2 Standard 2/3 Standard

1/4 Standard
1/2 Standard
2/3 Standard

**The next generation
of anterior matrices:**
We're taking the
guesswork out of
matrix selection



Contact your Bioclear Sales Rep to Pre-Order
WWW.BIOCLEARMATRIX.COM
1 800 770 5277





TSS

TOOTH & SURFACE SPECIFIC
MATRIX SYSTEM
BY BIOCLEAR



100 Anterior Matrices
100 Anterior & Premolar & Canine
Anterior & Premolar & Canine Matrices

75 Workups
100 Workup Matrices
100 Posterior Matrices

**The next generation
of anterior matrices:**
We're taking the
guesswork out of
matrix selection



Contact your BioClear Sales Rep to Pre-Order
www.bioclearmatrix.com
1 855 752 5327



TSS

TOOTH & SURFACE SPECIFIC
MATRIX SYSTEM
BY BIOCLEAR




100 Anterior Matrices
100 Wedges
100 Posterior Matrices
100 Posterior Wedges

**The next generation
of anterior matrices:**
We're taking the
guesswork out of
matrix selection.



Contact your Bioclear Sales Rep to Pre-Order
www.bioclearmatrix.com
1 855 712 5327




































Pre-operative



Immediate Post-Operative



One Week Post-Operative












Monolithic Kuraray Majesty
ES Paste & Flow: Shade
EW (extra white) with
Bioclear TSS Matrix









Day of treatment:
About 4 hours




1-week
post-op

Patient Chart

Patient Name



Patient ID




Date	Description	Provider	Tth	Surf	Status	Fee
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	6		Proposed	\$1,840.00
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	7		Proposed	\$1,840.00
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	8		Proposed	\$1,840.00
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	9		Proposed	\$1,840.00
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	10		Proposed	\$1,840.00
1/15/2026	2993A - BIOCLEAR 360 REJUVENATION COMPLEX	David J. Clark, D.D.S.	11		Proposed	\$1,840.00
1/15/2026	199.2 - BIOCLEAR ORTHO CORRECTION PER TOOTH	David J. Clark, D.D.S.	10		Proposed	\$306.00
1/15/2026	299.5 - BIOCLEAR DIASTEMA CLOSURE	David J. Clark, D.D.S.	10	MD	Proposed	\$285.00




Total Bioclear Case Fee?

\$11,631



The Dahl Technique



The Dahl "Light" Technique



The Epidemic of Cracked and Fracturing Teeth



Dr. David J. Frazier, D.D.S., M.S.D.

Tooth and denture failure is a real number. Cracked and fractured teeth are one of the most leading causes of tooth loss in industrialized nations.¹ Fractured teeth with conservative analgesic and antibiotic use, softening, and long-term crown placement, both anterior and posterior, are fracturing at all three levels. Fractured teeth treated with softening, softening, and softening (Figure 1). Additionally, endodontics are reporting that cracked teeth are now replacing crown teeth for sensitive etiologies in a number of patients referred for endodontic treatment.

We assume that the reason of the prevalence and growing problem is not so simple. Larger (3) our stress levels are rising, and chewing and tension are also more general, it is the advent of the high-speed handpieces, speeded billions of horsepower (10,000) every preparation that have work and powder teeth (Figure 2) and the deep cast crown preparations in the case of powder conductive, and yet aggressive "cross-link" radiolysis change in denture, and more teeth.

Dentistry has great options (including a poorly understood process). A dramatic change in the form of combining dental materials, dental, dentistry, blending, and implants has shifted our attention more from some of the most original aspects of the building and of dentistry, and dentistry, and dentistry. We now find our profession faced with a tremendous problem that will require a significant commitment and reevaluation of



Figure 1. A tooth with a crack and a filling. The tooth is cracked and the filling is cracked. The tooth is cracked and the filling is cracked.

resources to address appropriately. There is no need to assess. Many, to most simply will require a new part to work.

BRINGING DIAGNOSIS AND PREVENTION OFF OF THE SCENE

Microscopic and optical-drive diagnosis have been the accepted standards for the diagnosis of cracked teeth. The following list

of items of the lack of visual endodontics results therapies that education can be used in the treatment process. They feature first experience of vision through the clear advantage in the staggering array of cracks that can occur within such structures (Figure 3). Traditional visual exam (visual or biopsy) limits the dentist's ability to assess the pro-

cessed on page 48

Phase 1:	Phase 2:	Phase 3:	Phase 4:	Phase 5:
Diagnose the root of the problem (to be used in the treatment plan)	Reduce or eliminate the cause of the problem (to be used in the treatment plan)	Correct the problem (to be used in the treatment plan)	Prevent the problem (to be used in the treatment plan)	Re-evaluate the problem (to be used in the treatment plan)

1. Most teeth in aging adults display enamel cracks.
2. Cracked teeth, even traumatic cracks, do not necessarily indicate that the tooth is cracked.
3. Only enamel cracks do not penetrate significantly (70-80%).
4. Only enamel cracks that do not penetrate significantly (70-80%) are considered cracked.
5. Three types of cracking patterns are possible: enamel cracks, dentin cracks, and intermediate enamel (not contributing to microleakage enamel is radiolysis).
6. Dentin cracks should be considered as enamel cracks.
7. Dentin cracks that penetrate into the pulp (to be used in the treatment plan) are considered as enamel cracks (to be used in the treatment plan).
8. Only teeth with both types of enamel cracks (enamel cracks and dentin cracks) should be considered as structurally unsound.

Reprinted/Modified: Frazier, D.D.S., M.S.D.

“Dentistry’s last great mystery, that of fracturing, is a poorly understood process.”




**ONE OF THE ENDURING
FIRST IMPRESSIONS OF
CLINICAL MICROSCOPE
DENTISTRY IS THE
STAGGERING ARRAY OF
CRACKS THAT ARE
SUDDENLY VISIBLE**


Clark, D.J., Paquette J., Sheets C., *Journal of Esthetic and Restorative
Dentistry*; Special Edition, 2004

Nature does
nothing
uselessly.

Aristotle 322 BC




Enamel and dentin have radically different physical characteristics, as do their corresponding cracks.

A close-up photograph of a tooth, likely a maxillary premolar, showing a vertical crack that extends through the enamel and into the dentin. The crack is clearly visible as a dark line. The tooth is set against a dark background, and the lighting highlights the texture and color of the enamel and dentin. The text is overlaid on the upper portion of the image.


Enamel and dentin have radically different physical characteristics, as do their corresponding cracks.

12X



12X

Which enamel crack is "native", which crack is indicative of underlying pathology?



JOURNAL OF ESTHETIC AND RESTORATIVE DENTISTRY

Official Publication of the American Academy of Esthetic Dentistry,
Japan Academy of Esthetic Dentistry, International Federation of Esthetic Dentistry,
American Academy of Cosmetic and Adhesive Dentistry,
British Academy of Esthetic Dentistry, Dutch Academy of Esthetic Dentistry,
and the Scandinavian Academy of Esthetic Dentistry

www.bcdedcker.com

Definitive Diagnosis of Early Enamel and Dentinal Cracks Based on Microscopic Evaluation

DAVID J. CLARK, DDS,
DORLAND G. SMITH, DDT,
JACQUELINE B. QUINN, DDS

ABSTRACT

The diagnosis of cracked tooth and incomplete coronal fracture have historically been symptomatic. The dental operating microscope at 40x magnification can fundamentally change a clinician's ability to diagnose such conditions.

Clinicians have been observing cracks under extreme magnification to be merely cosmetic. Patients are becoming dissatisfied that can lead to appropriate treatment options symptomatic or asymptomatic tooth structure repair. Generally, many cracks are not structural and can lead to microleakage and over-restoration. Methodologic microscopic examination, an understanding of crack progression, and an appreciation of the types of cracks will guide a clinician in making appropriate decisions.

Teeth can have structural cracks in various stages. To date, diagnosis and treatment are very subjective and range of crack descriptions.

CLINICAL SIGNIFICANCE

This article gives new guidelines for recognition, classification, and treatment of cracked teeth based on the routine use of 40x magnification. The significance of enamel cracks as they relate to coronal cracks is detailed.

(J Esthet Restor Dent 16:XXX-XXX, 2004)


Microscopic and symptomatic diagnoses have been the accepted modalities for cracked teeth. The inherent limitations of the lack of visual confirmation result in therapies that often do not cure the treatment process. One limiting factor, impression of vision through a clinical microscope is the emerging array of cracks that color visible tooth structures. Traditional visualization (unaided

or ocular) cannot reliably distinguish the presence or severity of the majority of these cracks (Figure 1).

At extreme magnification levels (40x and greater), the translucent nature of enamel yields a wealth of information, such as color changes within the enamel may indicate early decay, microleakage, and a lack of structural integrity of dentin and

enamel. Being able to see previously invisible decay can aid restorative decisions to more appropriate early treatment of compromised teeth before desiccating, traumatic, pulp involvement, and periodontal breakdown occur. The value of early diagnosis of the structural breakdown of teeth will become even more significant with our aging population coupled with increasing tooth retention in this population.

¹President, Academy of Microscopic Oral and Esthetic Dentistry
²Consultant, American Academy of Esthetic Dentistry, Newport Beach, CA; Clinical Professor of Restorative Dentistry, USC School of Dentistry, Los Angeles, CA, USA
³Consultant, American Academy of Esthetic Dentistry, Newport Beach, CA; Associate Professor, Restorative Dentistry, USC School of Dentistry, Los Angeles, CA, USA




Vertical Crack End Stage:
Complete Tooth Fracture
(split tooth)

Endstage Oblique Fracture Cuspal Fracture





Enamel and dentin have radically different physical characteristics, as do their corresponding cracks.



12X

Which enamel crack is "native", which crack is indicative of underlying pathology?



A microscopic image of an enamel crack, showing a dark, irregularly shaped crack running through the enamel structure. The crack is surrounded by a lighter, more crystalline material. The overall image has a blueish tint.

Nomenclature and Classification System for Enamel Cracks

Type I: Have little or no risk of
underlying pathology


Type I: Have little or no risk of underlying pathology

- A) Craze lines
- B) Small vertical cracks
- C) Cracks that follow natural anatomic grooves
- D) Cracks with superficial stain penetration

Type I: Have little or no risk of underlying pathology

- F) Cracks Resulting from Polymerization Shrinkage of Composite
- G) Enamel Crackling in Aging, Thin, Cervical Enamel


Craze Lines



- Linear
- Perfectly vertical
- Never extend into dentin
- Don't house debris or have halos
- Can look dramatic but are never structural

Nomenclature and
classification of
enamel cracks


Craze Lines



Do not widen as
they extend
toward the
occlusal

Nomenclature and
classification of
enamel cracks

C. Cracks that follow anatomic grooves




Nomenclature and
classification of
enamel cracks

C. Cracks that follow natural anatomic grooves

D Superficially Stained Cracks


Nomenclature and
classification of
enamel cracks



Type II: Have moderate risk of underlying pathology

- A) V-shaped enamel ditching with no adjoining restoration often associated with a wear facet centered over an otherwise benign crack
- B) V-shaped enamel ditching with a prior restoration often associated with a wear facet centered over an otherwise benign crack
- C) Cracks that detour from or do not follow natural anatomic grooves

B. V shaped ditching



Nomenclature and
classification of
enamel cracks

C. Non linear cracks that detour from or do not follow natural anatomic grooves



Type II: Have moderate risk of underlying pathology


- A) V-shaped enamel ditching with or without an adjoining restoration often associated with a wear facet centered over an otherwise benign crack
- B) Cracks that detour from or do not follow natural anatomic grooves




8X

C. Non linear cracks that
detour from or do not follow
natural anatomic grooves

12X






- Unusual or asymmetric gapping of filling material






C. Nonlinear cracks that detour from or do not follow natural anatomic grooves








Exercise 1
24X









Type III: Have high risk of underlying pathology

- A) Diagonal cracks branching off a vertical crack
- B) Horizontal or diagonal cracks that emanate from the corner of a restoration
- C) Cracks that house debris
- **D) The tooth looks gray**



A. Diagonal cracks



B. Diagonal cracks that emanate from the corner of a restoration




Nomenclature and
classification of
enamel cracks




A Two Year Study 2001-2003

- 120 Teeth
- 24 months
- Every tooth treated restoratively examined at 16x after restorative material removed




A Two Year Study

- 120 Teeth
- 188 Dentin Cracks
(Incomplete Fractures)



A Two Year Study



- The third leading cause of tooth loss today is from splits/fractures. (Lynch CD McCinnel RJ. The cracked tooth Syndrome. J Dent Assoc Sept 2002; 68(8):470-475)
- Crack initiation and progression should not be labeled a “syndrome”
- Diagnosis and treatment should be similar to approaches we utilize to address other pathologies (i.e. caries, periodontal disease)

B. Diagonal cracks that emanate from the corner of a restoration



2X

Nomenclature and
classification of
enamel cracks

B. Diagonal cracks that emanate from the corner of a restoration



4X

B. Diagonal cracks that emanate from the corner of a restoration



8X

B. Diagonal cracks that emanate from the corner of a restoration



B. Diagonal cracks that emanate from the corner of a restoration



12X





16X

B. Horizontal cracks



C. Crack that houses debris






Type III: Have high risk of underlying pathology

- D. Pair of cracks that outline an area (cusp or marginal ridge) of discolored enamel
- E. Crack with corresponding halo of brown, grey or white centered on crack

D. Pair of cracks that outline an area (cusp or marginal ridge) of discolored enamel



E. Crack with corresponding halo of brown, grey or white centered on crack



E. Crack with corresponding halo of brown, grey or white centered on crack



Other Microscopic Findings That Can Indicate a Lack of Coronal Structural Integrity

- Unusual or asymmetric gapping of filling material
- Crack in Filling material
- Gray discoloration of a cusp or cusps






A Tale of Three Marginal Ridges



A Tale of Three Marginal Ridges









A Tale of Three Marginal Ridges



A Tale of Three Marginal Ridges








A Tale of Three Marginal Ridges





A Tale of Three Marginal Ridges







What
we
expected




















Incomplete and Complete Fractures



Case Study

Occlusal view of distal






Initial
appointment



3-6 month
wait with
temp crown



4 month
wait





The Old Cracked Tooth Classification System...

Was Based on Symptoms
and Conjecture

Epidemiologic
third most com
countries. This
syndrome is of

Cracked tooth syndrome

Figure 1






Figure 1



- Sad
main



cold food or drink. The problem with waiting for symptoms to develop before a diagnosis is made is that at this point the crack is usually at what is referred to by Clark⁵ as end-stage crack progression, requiring more aggressive treatment.

The emergence of the dental microscope has been a tremendous aid in the diagnosis and treatment of cracked teeth. The routine, methodical

“Build a bridge
from buccal to
lingual, and the
crack(s)
becomes
dormant.”



DR. ALEX FOK



DC BIOCLEAR LEARNING CENTER

Tacoma WA. USA · Solihull UK

Varberg Sweden · Cairo Egypt

Syracuse Italy · Taubate Brazil

Livermore CA (Bioclear pediatrics)

Seoul Korea · Madrid/Barcelona

Sydney Australia · Baghdad Iraq

Getting paid to be conservative

BIOCLEAR VS. CROWNS AND VENEERS

Bioclear is an alternative to traditional methods for enhancing a smile. Rather than preparing for a crown or veneer, the Bioclear Method preserves the natural tooth structure, enamel, and tooth durability.




AN HONEST LOOK AT CROWN PREPARATIONS

As illustrated below, crowns, veneers, and onlays require the removal of a significant amount of healthy tooth structure. Bioclear dentists can leave most or all of the tooth structure. Bioclear is a very attractive option to patients.



~3/4 Crown
47% of tooth removed

Full Crown
76% of tooth removed



BIOCLEAR PREPARATION

Bioclear allows dentists to **conserve** healthy tooth structure



CROWN PREPARATION

Crowns require dentists to **remove an average of 76%** of the tooth structure prior to the procedure



VENEER PREPARATION

Veneers require dentists to **remove an average of 47%** of the tooth structure prior to the procedure

∞ BIOCLEAR

Stop by convention booth #14 to...

- ✓ Hear more about Bioclear courses
- ✓ Meet our team of Bioclear experts
- ✓ Order products
- ✓ Register for courses
- ✓ Learn why Bioclear is a BIG DEAL!




Treat yourself to a weekend with  BIOCLEAR

BLACK TRIANGLE

CERTIFICATION COURSE

18 CE CREDITS

This three-part live hands-on certification course will teach you how to treat black triangles, gingival recession, root abrasions, and perform confident restorations. In becoming a certified Bioclear black triangle doctor, you'll increase your overall skill and knowledge of Bioclear and learn to market your new skills to patients.



Upcoming Dates:

Denver

October 9th 2026

Walnut Creek

August 21st 2026

Limited spots available!

Scan to learn more
& sign up for your
local BT Course!



BEFORE



AFTER



BEFORE



AFTER



∞ BIOCLEAR LEARNING CENTER

The Bioclear Learning Center aims to improve restorative outcomes and raise patient expectations of modern composite dental care by providing dentists with continuing education and certification in the practice of the Bioclear Method.

- The Learning Center is equipped with exclusive, patented Bioclear tools, equipment, and your own operator setup
- Small class sizes with hands-on, interactive curriculum modules
- Online class options available at the convenience of your practice
- Operate a Tell, Show, Do education model
- Learn and understand the Five Pillars of the Bioclear Method

4-DAY CORE ANTERIOR & POSTERIOR COURSE

36 CE CREDITS

This course builds a strong foundation in the Bioclear Method and an understanding of Bioclear products. The focus of the course is indirect methods on adult dentitions as an alternative procedure to porcelain crowns and veneers.

Students will learn the foundation of posterior and anterior restorations and are introduced to the engineering principles involved in the design of the new non-retentive compression-based preps. Students collaboratively practice all applications of the Bioclear Method during intensive, hands-on exercises that simulate posterior and anterior restorations.



LEARN MORE



BEFORE






AFTER















2026 Certification Courses



- Core Anterior + Core Posterior • Complex Cases & Problem Solving
- Smile Design & Comprehensive Anterior Rejuvenation • Bioclear Alumni Summit

January

S	M	T	W	T	F	S
				1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31

February

S	M	T	W	T	F	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28

March

S	M	T	W	T	F	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

April

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30		

May

S	M	T	W	T	F	S
					1	2
3	4	5	6	7	8	9
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

June

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30				

July

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

August

S	M	T	W	T	F	S
						1
2	3	4	5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31					

September

S	M	T	W	T	F	S
		1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30			

October

S	M	T	W	T	F	S
				1	2	3
4	5	6	7	8	9	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31

November

S	M	T	W	T	F	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30					

December

S	M	T	W	T	F	S
		1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30	31		

2027 Certification Courses



- Core Anterior + Core Posterior • Complex Cases & Problem Solving
- Smile Design & Comprehensive Anterior Rejuvenation

JANUARY

S	M	T	W	T	F	S
					01	02
03	04	05	06	07	08	09
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

FEBRUARY

S	M	T	W	T	F	S
	01	02	03	04	05	06
07	08	09	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28						

MARCH

S	M	T	W	T	F	S
	01	02	03	04	05	06
07	08	09	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

APRIL

S	M	T	W	T	F	S
				01	02	03
04	05	06	07	08	09	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	

MAY

S	M	T	W	T	F	S
						01
02	03	04	05	06	07	08
09	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31					

JUNE

S	M	T	W	T	F	S
		01	02	03	04	05
06	07	08	09	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30			

JULY

S	M	T	W	T	F	S
				01	02	03
04	05	06	07	08	09	10
11	12	13	14	15	16	17
18	19	20	21	22	23	24
25	26	27	28	29	30	31

AUGUST

S	M	T	W	T	F	S
01	02	03	04	05	06	07
08	09	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

SEPTEMBER

S	M	T	W	T	F	S
			01	02	03	04
05	06	07	08	09	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30		

OCTOBER


S	M	T	W	T	F	S
					01	02
03	04	05	06	07	08	09
10	11	12	13	14	15	16
17	18	19	20	21	22	23
24	25	26	27	28	29	30
31						

NOVEMBER


S	M	T	W	T	F	S
	01	02	03	04	05	06
07	08	09	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30				



DECEMBER


S	M	T	W	T	F	S
				01	02	03
				04	05	06
07	08	09	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			





Dr. Rafael Bustamante
Full Faculty at the
Bioclear Learning center







- Selective Caries Removal
- Calla Lilly Compression Based Cavity Preparation for a Cracked Tooth
- Deep Margin Acquisition
- Dealing with a Furcation




This is *Dr. Bustamante's* finger
BTW






What changed between pre and post rubber dam photos?





















The new “Long Diamond Wedge” Available
in all five Diamond Wedge shapes



Long (Double-Hump) Furcal Wedge







Outline of Radius Bevel

Before Blasting





After Blasting and Final
Radius Beveling




















Miles ahead of everybody else for Bulk Fills






YIKES!!


GC “Bulk Injectable”

Discussion:
Polymerization
shrinkage stress and
cuspal deflection is
very real but
challenging to
explain




Discussion:

The volume of composite, and residual strength of the tooth are key players. If half of tooth is missing, ie there is a good sized hole in the occlusal or interproximal, we are recommending Filtek One (bulk fill paste) plus GC bulk flow or Filtek bulk flow.



Discussion:

If you we have a fair amount of dentin removed, we need Filtek One



Miles ahead of everybody else for Bulk Fills







DC BIOCLEAR LEARNING CENTER

Tacoma USA · Solihull UK
Varberg Sweden · Cairo Egypt
Syracuse Italy · Taubate Brazil
Livermore CA (Bioclear pediatrics)
Seoul Korea · Madrid/Barcelona
Sydney Australia · Baghdad Iraq








Last week's attendees at the certification course








Hands-On Workshop: Better, Faster, Prettier Composite Dentistry - 3 CE 1:45PM-4:15PM



- Audience: Hands-On Workshop Pre-registered \$100
- Location: Room 410 C, Boise Centre East
- **by Dr. David Clark, Bioclear**

Experience a fresh approach to modern composite dentistry with Dr. Clark. Say goodbye to traditional techniques like layering, condensing, and burnishing. Join us to master injection molding, achieving precise margins, smooth contours, and flawless, mirror-like restorations in both black triangle and Class II scenarios.

- Learn: The 6 steps to achieve rock solid posterior contacts in even the most difficult cases
- Learn: How can you combine flowable & regular composites with the injection molding technique & when to use bulk fills
- Learn: Treating the dreaded black triangle, especially common after tray aligner therapy
- Learn: How to handle heated composites

 solventum
(3M)

3M™ Filtek™ One
Bulk Fill Restorative


Scientific
presentation



Innovative methacrylate monomers for lower shrinkage and proven stress relief

AUDMA: Aromatic urethane dimethacrylate

- Higher molecular weight with less number of reactive groups
- Moderates volumetric shrinkage
- Contributes to stress reduction




AFM: Addition-fragmentation (AF) monomer

- Reacts into developing polymer network through terminal methacrylate bonds like other dimethacrylate monomers
- Central AF group can fragment and release stress
- Fragment may then polymerize into network in a lower stress orientation compared to its pre-fragmented state



Polymerization shrinkage


Comparison of common bulk fill composites



Source: Solventum Internal Data

3M™ Filtek™ One Bulk Fill Restorative


is equivalent to 3M™ Filtek™ Bulk Fill Posterior and has lower polymerization shrinkage than Venus® Bulk Fill, Surefil® SDR® Flow +.



Bonded disc method


Polymerization stress

Comparison of common universal composites



3M™ Filtek™ One Bulk Fill Restorative


has similar polymerization stress when placed in bulk compared to many of the universal materials that are placed in increments.



Source: Solventum Internal Data


Polymerization stress

Comparison of common bulk fill composites



3M™ Filtek™ One Bulk Fill Restorative

has similar polymerization stress compared to 3M™ Filtek™ Bulk Fill Posterior Restorative and SonicFill™ 2.



Source: Solventum Internal Data